

## Futura - 1 Spill



Ardent Futura is a high copper non gamma-2 amalgam alloy based on a mixture of irregular atomized and lathe cut particles, especially designed for outstanding corrosion resistance and lowest creep value. Technical Data: In accordance to ISO 24234 and ISO 13897 • Compressive strength 1h: 210 MPa. • Compressive strength 24h: 530 MPa. • Static creep: 0,05%. • Dimensional change: +0,02%. • Working time: ~4 min. Delivery forms • Standard package of 50 capsules in one, two or three spill.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)